




**Sterilization, Miscellaneous**  
**Sterilisation, Verschiedenes**  
**Esterilización, Diversos artículos**  
**Stérilisation, Divers articles**  
**Sterilizzazione, Articoli**



DAVIS

43-100-16 16 cm, 6¼"  
43-100-24 24 cm, 9½"



43-104-16  
16,5 cm, 6½"



UNIVERSAL

43-110-20 20 cm, 8" 43-111-20 20 cm, 8"  
43-110-28 28 cm, 11" 43-111-28 28 cm, 11"

Curved  
Kappe  
curva  
courbe  
curva



ROGGE  
43-114-20  
20 cm, 8"

ROGGE  
43-114-30  
30 cm, 12"




ROGGE  
43-118-23  
23 cm, 9"



CHEATTLE  
43-124-27  
27 cm, 10½"



CHEATTLE  
43-126-28  
28 cm, 11"



**43-130-00**

Sterilizing forceps consisting of:

- 43-130-01 1 Sterilizing forceps 25 cm, curved, with lid
- 43-130-02 1 Jar with wide bottom plate (anti-tilt)

Kappenzange bestehend aus:

- 43-130-01 1 Kornzange 25 cm, gebogen, mit Kappe (Deckel)
- 43-130-02 1 Standgefäß mit großer Bodenplatte (standsicher)

Combinación tapa/pinza compuesta de:

- 43-130-01 1 Pinza para vendaje 25 cm, curva, con tapa
- 43-130-02 1 Tarro con placa de pie grande (estable)

Combinaison couvercle/pince, composée de:

- 43-130-01 1 Pince à tampons 25 cm, courbe, avec couvercle
- 43-130-02 1 Eprouvette porte-pince avec large pied (estable)

Sterilizzatore per pinze composto di:

- 43-130-01 1 Pinza per medicatura 25 cm, curva, con coperchio
- 43-130-02 1 Contenitore con fondo largo (stabile)



**43-132-00**

Sterilizing forceps consisting of:

- 43-132-01 1 Sterilizing forceps 32 cm, curved, with lid
- 43-132-02 1 Jar with wide bottom plate (anti-tilt)

Kappenzange bestehend aus:

- 43-132-01 1 Kornzange 32 cm, gebogen, mit Kappe (Deckel)
- 43-132-02 1 Standgefäß mit großer Bodenplatte (standsicher)

Combinación tapa/pinza compuesta de:

- 43-132-01 1 Pinza para vendaje 32 cm, curva, con tapa
- 43-132-02 1 Tarro con placa de pie grande (estable)

Combinaison couvercle/pince, composée de:

- 43-130-01 1 Pince à tampons 25 cm, courbe, avec couvercle
- 43-130-02 1 Eprouvette porte-pince avec large pied (estable)

Sterilizzatore per pinze composto di:

- 43-132-01 1 Pinza per medicatura 32 cm, curva, con coperchio
- 43-132-02 1 Contenitore con fondo largo (stabile)




**43-136-58**  
58 cm, 22 3/4"



**SAUERBRUCH**  
**43-138-55**  
55 cm, 21 1/2"

**Pick-up forceps**  
Abfallzange  
Pinza para desechos  
Pince à déchets  
Pinza raccogliatrice



**BUNT**  
 43-140-10  
 10,5 cm, 4 1/4"

Safety pin  
 Sterilisierklammer  
 Imperdible para esterilización  
 Epingle à instruments  
 Spillo di sicurezza



Safety pin  
 Sterilisierklammer  
 Imperdible para esterilización  
 Epingle à instruments  
 Spillo di sicurezza

**MAYO**  
 43-142-14  
 14cm, 5 1/2"




Safety pin  
 Sterilisierklammer  
 Imperdible para esterilización  
 Epingle à instruments  
 Spillo di sicurezza

**MAYO**  
 43-144-14  
 14 cm, 5 1/2"

**Safety rack for 36 ring handle instruments, 1 Pair**  
 Sterilisiergestell für 36 Ringgriffinstrumente, 1 Paar  
 Porta-instrumento para 36 instrumentos con mangos de anillo, 1 par  
 Porte instruments pour 36 instruments à anneaux, 1 paire  
 Telaio per 36 strumenti, 1 paio

**43-146-15**  
 15,5 cm, 6"



**03-146-25**

25,5 cm, 10"  
**Safety rack for 60 ring handle instruments, 1 Pair**  
 Sterilisiergestell für 60 Ringgriffinstrumente, 1 Paar  
 Porta-instrumento para 60 instrumentos con mangos de anillo, 1 par  
 Porte instruments pour 60 instruments à anneaux, 1 paire  
 Telaio per 60 strumenti, 1 paio



a = Gesamtänge  
 longueur totale  
 lunghezza totale  
 b = Bürstenlänge  
 longueur de brosse  
 lunghezza delle spazzole  
 c = Bürstendurchmesser  
 diamètre de brosse  
 diametro delle spazzole

43-147-03

43-149-03

43-151-07

43-151-08

43-151-13

43-152-13

a Abmessungen in mm  
 b dimensions en mm  
 c dimensioni in mm.



43-153-08



43-153-35



43-153-13



43-153-25



43-153-28

A	215 mm	215 mm	215 mm	215 mm	215 mm
B	60 mm	65 mm	95 mm	95 mm	100 mm
C	13 mm	15 mm	20 mm	25 mm	28 mm



43-155-12



43-155-50



43-155-35



43-155-60



43-157-50  
330 x 130 x 50 mm



43-157-95  
210 x 95 mm



KONRICH



SCHIMMELBUSCH

Ø cm	Height Höhe Altura hauteur altezza		Height Höhe Altura hauteur altezza		Filters only Filter, allein Filtros sueltos Filtres seuls Filtri, solo		
	▼ cm	Ø in	▼ in	KONRICH	KONRICH	SCHIMME- LBUSCH	
8	8	3¼	3¼	43-160-01	43-161-01	43-166-01	
10	10	4	4	43-160-02	43-161-02	43-166-02	
12	12	4¼	4¼	43-160-03	43-161-03	43-166-03	
14	14	5½	5½	43-160-04	43-161-04	43-166-04	
16	12	6¼	4¼	43-160-05	43-161-05	43-166-05	
16	15	6¼	6	43-160-06	43-161-06	43-166-06	
19	15	7½	6	43-160-07	43-161-07	43-166-07	
19	19	7½	7½	43-160-08	43-161-08	43-166-08	
24	12	9½	4¼	43-160-09	43-161-09	43-166-09	
24	16	9½	6¼	43-160-10	43-161-10	43-166-10	
24	24	9½	9½	43-160-11	43-161-11	43-166-11	
29	16	11½	6¼	43-160-12	43-161-12	43-166-12	
29	24	11½	9½	43-160-13	43-161-13	43-166-13	
34	16	13½	6¼	43-160-14	43-161-14	43-166-14	
34	24	13½	9½	43-160-15	43-161-15	43-166-15	
39	14	15½	5½	43-160-16	43-161-16	43-166-16	
39	19	15½	7½	43-160-17	43-161-17	43-166-17	
39	29	15½	11½	43-160-18	43-161-18	43-166-18	



- 43-170-01 260 x 250 x 50 mm
- 43-170-02 260 x 250 x 75 mm
- 43-170-03 260 x 250 x 100 mm

**Sterilizing trays of perforated steel sheet, for container 300 mm long**

Sterilisier-Siebschalen aus perforiertem Blech, für Container 300 mm lang

Bandejas de esterilización perforadas de chapa de acero para contenedores de 300 mm de largo

Plateaux en inox perforé pour stérilisation, pour container de 300 mm de longueur

Castelli in acciaio perforato per la sterilizzazione di contenitori lunghi 300 mm



- 43-174-01 480 x 250 x 50 mm
- 43-174-02 520 x 250 x 50 mm
- 43-174-03 520 x 250 x 80 mm
- 43-174-04 520 x 250 x 100 mm


**Sterilizing trays of perforated steel sheet, for container 600 mm long**

Sterilisier-Siebschalen aus perforiertem Blech, für Container 600 mm lang

Bandejas de esterilización perforadas de chapa de acero para contenedores de 600 mm de largo

Plateaux en inox perforé pour stérilisation, pour container de 600 mm de longueur


Castelli in acciaio perforato per la sterilizzazione di contenitori lunghi 600 mm



¼

- 43-178-01 215 x 160 x 18 mm
- 43-178-02 235 x 180 x 15 mm
- 43-178-03 300 x 150 x 20 mm
- 43-178-04 385 x 210 x 10 mm
- 43-178-05 350 x 240 x 10 mm
- 43-178-06 400 x 270 x 20 mm
- 43-178-07 420 x 280 x 10 mm


**Instrument storing trays**  
Instrumentenablageplatten  
Bandejas para instrumentos  
Plateaux pour instruments  
Vassoi porta strumenti



¼


- 43-180-01 230 x 155 x 30 mm
- 43-180-02 300 x 175 x 30 mm
- 43-180-03 300 x 175 x 40 mm
- 43-180-04 350 x 240 x 36 mm
- 43-180-05 400 x 300 x 50 mm

**Universal trays**  
Mehrzweckschalen  
Bandejas universales  
Bassins universals  
Bacinele universali




**43-184-01**  
235 x 190 x 40 mm

**Tray**  
Schwenschale  
Bandeja rectangular  
Bassin et cuvette pour guéridons  
Bacinelle rettangolari




**Lid without knob**  
Deckel ohne Knopf  
Tapa sin botón  
Avec couvercle sans bouton  
Con coperchio senza bottone

<b>43-186-01</b>	165 x 85 x 35 mm
<b>43-186-02</b>	200 x 100 x 35 mm
<b>43-186-03</b>	200 x 100 x 50 mm
<b>43-186-04</b>	230 x 130 x 50 mm
<b>43-186-05</b>	260 x 150 x 50 mm
<b>43-186-06</b>	270 x 235 x 50 mm
<b>43-186-07</b>	300 x 200 x 50 mm
<b>43-186-08</b>	325 x 275 x 50 mm



**Lid with knob**  
Mit Knopfdeckel  
Tapa con botón  
Avec couvercle à bouton  
Con coperchio a bottone

<b>43-186-01</b>	165 x 85 x 35 mm
<b>43-180-02</b>	200 x 100 x 35 mm
<b>43-180-03</b>	200 x 100 x 50 mm
<b>43-180-04</b>	220 x 130 x 50 mm
<b>43-180-05</b>	260 x 150 x 50 mm
<b>43-180-06</b>	270 x 235 x 50 mm
<b>43-180-07</b>	300 x 200 x 50 mm
<b>43-180-08</b>	325 x 275 x 50 mm



**43-190-01**  
205 x 105 x 35 mm

**Lid without knob**  
Deckel ohne Knopf  
Tapa sin botón  
Avec couvercle sans bouton  
Con coperchio senza bottone



**Round bowls**  
Runde Schalen  
Cubetas redondas  
Cupules  
Scodelle rotonde

		<b>Capacity</b> Inhalt Cabida Capacité Capacità
<b>43-200-01</b>	40 x 19 mm	0,05 L
<b>43-200-02</b>	61 x 30 mm	0,07 L
<b>43-200-03</b>	80 x 40 mm	0,14 L
<b>43-200-04</b>	116 x 50 mm	0,35 L
<b>43-200-05</b>	128 x 55 mm	1,45 L
<b>43-200-06</b>	147 x 65 mm	1,70 L
<b>43-200-07</b>	157 x 75 mm	1,10 L
<b>43-200-08</b>	187 x 85 mm	1,60 L



**Wash bowls**  
Waschschüsseln  
Palanganas  
Cuvettes  
Scodella per liquidi

	<b>Sizes (Ø x H)</b> Abmessungen (Ø X H) Tamaños (Ø X a) Dimensions (Ø X la) Dimensioni (Ø X h)	<b>Capacity</b> Inhalt Cabida Capacité Capacità
<b>43-204-01</b>	300 x 90 mm	6,0 L
<b>43-204-02</b>	330 x 115 mm	0,0 L



¼

**Solution bowls**  
Spülschalen  
Cubetas para soluciones  
Cuvettes  
Scodelle per soluzioni

**Sizes (∅ x H)**  
Abmessungen (∅ X H)  
Tamaños (∅ X a)  
Dimensions (∅ X la)  
Dimensioni (∅ X h)

**Capacity**  
Inhalt  
Cubida  
Capacité  
Capacità

43-240-01	70 x 35 mm	0,1 L
43-240-02	90 x 50 mm	0,2 L
43-240-03	100 x 65 mm	0,3 L
43-240-01	115 x 65 mm	0,4 L
43-240-02	120 x 70 mm	0,5 L
43-240-03	155 x 85 mm	1,7 L
43-240-04	195 x 120 mm	2,0 L
43-240-05	220 x 120 mm	2,7 L
43-240-06	245 x 130 mm	4,0 L
43-240-07	290 x 150 mm	5,5 L
43-240-08	305 x 175 mm	8,0 L
43-240-09	360 x 185 mm	12,0 L
43-240-10	370 x 200 mm	14,0 L




½



½

**Sizes (∅ x H)**  
Abmessungen (∅ X H)  
Tamaños (∅ X a)  
Dimensions (∅ X la)  
Dimensioni (∅ X h)

43-250-01	30 x 90 mm	43-252-01
43-250-02	50 x 100 mm	43-252-02
43-250-03	50 x 130 mm	43-252-03
43-250-04	50 x 175 mm	43-252-04
43-250-05	75 x 100 mm	43-252-05
43-250-06	75 x 130 mm	43-252-06
43-250-07	75 x 175 mm	43-252-07
43-250-08	100 x 175 mm	43-252-08
43-250-09	100 x 200 mm	43-252-09



¾



¾


43-260-01	75 x 60 mm
43-260-02	75 x 75 mm
43-260-03	100 x 60 mm
43-260-04	100 x 100 mm
43-260-05	125 x 80 mm
43-260-06	125 x 125 mm
43-260-07	150 x 150 mm

**Containers, lid with knob**  
Dosen, mit knopfdeckel  
Recipientes, tapa con botón  
Récipients avec couvercle à bouton  
Recipienti con coperchio a bottone




43-262-15






1/4

Kidney bowls	Sizes (l x w x h)
Nierenschalen	Abmessungen (l x b x h)
Cuñas	Tamaños (l x l x a)
Bassins réniformes	Dimensions (lo. x hau. x la.)
Bacinelle reniformi	Dimensioni (L x l x h)
43-300-17	170 X 100 X 35 mm
43-300-20	200 X 100 X 40 mm
43-300-22	250 X 140 X 40 mm
43-300-27	275 X 150 X 45 mm




1/4

43-310-23 Ø 235 mm  
43-310-31 Ø 310 mm  
**Bed pan ,lid with knob**  
Stechbecken, mit Knopfdeckel  
Cubetas, tapa con botón  
Bassin de lit  
Padelle, coperchio con bottone



1/4

43-316-23  
350 x 290



1/4

43-320-01  
**Without lid**  
ohne deckel  
sin tapa  
sans couvercle  
senza coperchio

43-320-02  
**With lid**  
mit deckel  
con tapa  
avec couvercle  
con coperchio



1/4

43-330-01  
**Triangular basin**  
Dreieckschale  
Cubeta triangulares  
Bassin triangulaires  
Bacinelle triangolari



1/4

43-330-02  
**Triangular basin**  
Dreieckschale  
Cubeta triangulares  
Bassin triangulaires  
Bacinelle triangolari



1/2

43-350-20  
25 cc

**Medicine cup**  
Medizinbecher  
Vaso graduado  
Gobelet  
Bicchiere per medicinali



1/2

43-354-50  
50 cc

**Medicine cup**  
Medizinbecher  
Vaso graduado  
Gobelet  
Bicchiere per medicinali



1/2

43-358-25  
20 cc

**Medicine spoon**  
Medizinlöffel  
Cuchara graduada  
Cuillère pour médicaments  
Cucchiao per medicinali



1/2

43-360-01

**Bowl for dentures**  
Zahnprothessenschale  
Tarro para prótesis dental  
Réceptif à prothèses dentaires  
Recipienti per protesi dentarie



1/2

43-370-01  
43-370-02  
43-370-03

**Measuring jugs**  
Messbecher  
Copa graduada  
Gobelets gradués  
Misurino graduato



1/2

43-372-01  
43-372-02  
43-372-03  
43-372-04  
43-372-05

**Measuring jugs**  
Messbecher  
Copa graduada  
Gobelets gradués  
Misurino graduato



1/2

**Funnels**  
Trichter  
Embudos  
Entonnoirs  
Imbuti

43-380-01 100 mm  
43-380-02 120 mm  
43-380-03 150 mm  
43-380-04 200 mm  
43-380-05 240 mm



1/2

43-384-01

**Sputum cup**  
Sputumbecher  
Escupideras  
Crachoir  
Sputacchiera



1/2

43-388-01  
43-388-02  
43-388-03

**Irrigators**  
Irrigatoren  
Irrigadores  
Irrigateurs  
Irrigatori