









Tracheotomy
Tracheotomie
Traqueotomia
Trachéotomie
Tracheotomia








Fig.	0	1	2	3	4	5	6	7	8	9
Ø mm	④	⑤	⑥	⑦	⑧	⑨	⑩	⑪	⑫	⑬



LUER
41-300-00
41-300-09

Fig.	ø mm	L mm		Radius	Fig.	ø mm	L mm		Radius
0	4	45	41-300-00	N	5	9	70	41-300-06	N
1	5	50	41-300-01	N	6	10	75	41-300-07	N
2	6	55	41-300-02	N	7	11	80	41-300-08	N
3	7	60	41-300-03	N	8	12	85	41-300-09	N
4	8	65	41-300-04	N	9	13	90		



LUER
41-310-04
41-310-06

Fig.	ø mm	L mm		Radius
4	8	95	41-310-04	M
5	9	95	41-310-05	M
6	10	95	41-310-06	M
7	11	95	41-310-07	M
8	12	95	41-310-08	M
9	13	95	41-310-09	M



CHEVALIER
JACKSON
41-320-00
41-320-09

outer tube closed
inner tube closed

Außenrohr geschlossen
Innerrohr geschlossen

tubo exterior cerrado
tubo interior cerrado

canule externe

Fig.	ø mm	L mm		Radius
0	4	45	41-320-00	N
1	5	50	41-320-01	N
2	6	55	41-320-02	N
3	7	60	41-320-03	N
4	8	65	41-320-04	N
5	9	70	41-320-05	N
6	10	75	41-320-06	N
7	11	80	41-320-07	N
8	12	85	41-320-08	N
9	13	90	41-320-09	N



Regular radius
Normalradius
Radio normal
Rayon normal
Raggio normale

Fig.	ø mm	L mm	German silver Neusilber Plata meneses Argent blanc Argento bianco
4	8	65	41-330-04
5	9	70	41-330-05
6	10	75	41-330-06
7	11	80	41-330-07
8	12	85	41-330-08
9	13	90	41-330-09

Valve speaking tubes, with detachable round valve
Ventil-Sprechkanülen mit rundem aufsteckventil
Cánulas traqueales con válvula para hablar, redonda
Canules parlantes avec soupape ronde amovible
Cannule parlante con valvola rotonda mobile



Regular radius
Normalradius
Radio normal
Rayon normal
Raggio normale

Fig.	ø mm	L mm	Silver Silber Plata Argent Argento
4	8	65	41-340-04
5	9	70	41-340-05
6	10	75	41-340-06
7	11	80	41-340-07
8	12	85	41-340-08
9	13	90	41-340-09

Valve speaking tubes, inner tube with sliding valve

Ventil-Sprechkanülen, Innenrohr mit schiebeventil
Cánulas traqueales con válvula para hablar, cánula interior con válvula de corredera
Canules parlantes, tube intérieur avec soupape à coulisse
Canule parlanti, tubo interno con valvola scorrevole



Medium radius
Mittlerer Radius
Radio medio
Rayon moyen
Raggio medio

Fig.	ø mm	L mm	German silver Neusilber Plata meneses Argent blanc Argento bianco
4	8	95	41-350-04
5	9	95	41-350-05
6	10	95	41-350-06
7	11	95	41-350-07
8	12	95	41-350-08
9	13	95	41-350-09

Valve speaking tubes, inner tube with sliding valve

Ventil-Sprechkanülen, Innenrohr mit schiebeventil
Cánulas traqueales con válvula para hablar, cánula interior con válvula de corredera
Canules parlantes, tube intérieur avec soupape à coulisse
Canule parlanti, tubo interno con valvola scorrevole

ø 12 mm



41-360-26
26 cm 10 1/4"


Cleaning brush for trachea tubes

Reinigungsbürste für trachealkanülen

Cepilo de limpieza para cánulas traqueales

Brosse à nettoyer les canules trachéales


Spazzolino per pulire per cannule tracheali



41-370-00

Tracheotomy set, complete, consisting of:

- 1 Rack
- 1 Metal case
- 2 Handles No. 4
- 1 Blade fig. 20
- 1 Blade fig. 23
- 2 PEAN Artery forceps, 13 cm, 5"
- 1 Dissecting forceps, 13 cm, 5"
- 2 Tissue forceps, 13 cm, 5" 1x2 teeth
- 1 Grooved director, 13 cm, 5"
- 2 Trachea retractors, blunt, 16,5 cm, 6½"
- 1 FINSSEN retractor, blunt, 7 cm, 2¾"
- 1 BOSE trachea retractor, sharp, left, 16 cm, 6¼"
- 1 BOSE trachea retractor, sharp, right, 16 cm, 6¼"
- 1 LUER trachea tube, fig. 3, ø 7 mm
- 1 LUER trachea tube, fig. 4, ø 8 mm
- 1 Cleaning brush, 30 cm, 12"



41-380-00

Tracheotomy set, complete, consisting of:

- 41-380-01 1 Rack
- 41-380-02 1 Rack
- 41-380-03 1 Metal case
- 05-162-04 1 Handle No. 4
- 05-150-20 1 Blade fig. 20
- 05-150-22 1 Blade fig. 22
- 05-150-23 2 Blades fig. 23
- 11-184-14 4 CRILE Artery forceps, straight, 14 cm, 5½"
- 09-260-13 1 Dissecting forceps, 13 cm, 5"
- 09-461-13 2 Tissue forceps, 13 cm, 5½" 1x2 teeth
- 19-132-13 1 Grooved director, 13 cm, 5½"
- 17-216-01 2 Trachea retractors, blunt 16,5 cm, 6½"
- 41-120-16 1 BOSE Trachea retractor, sharp, left, 16,5 cm, 6½"
- 41-122-16 1 BOSE Trachea retractor, sharp, right, 16,5 cm, 6½"
- 41-200-14 1 TROUSSEAU Trachea dilator, 14 cm, 5½"
- 41-300-02 1 LUER Trachea tube, fig. 2, ø 6 mm, German silver
- 41-300-04 1 LUER Trachea tube, fig. 4, ø 8 mm, German silver
- 41-360-30 1 Cleaning brush, 30 cm, 12"